

FAX OR E-MAIL COMPLETED FORMS

Fax 866.700.7837 or e-mail to orders@surestep.net
Customer Service 877.462.0711

Date _____ Right Left Bilateral

Pt Last Name _____

Pt First Name _____

DOB _____ Weight _____ Sex M F

Dx _____

Company Name _____

Contact _____

Phone Number _____

Fax Number _____

Billing Address _____

City _____ ST _____ Zip _____

Shipping Address _____

City _____ ST _____ Zip _____

E-mail _____

PO _____

Order Confirmation: Fax E-mail

Device Style Options

Custom FO Right Left Bilateral

Custom UCBL Right Left Bilateral

Heel Post Dorsum Wings

Trimlines: Traditional Surestep

Elastic Strap

Custom UCBL with FO Right Left Bilateral

Heel Post Dorsum Wings

Trimlines: Traditional Surestep

Elastic Strap

Foot Orthotic Pattern Options

Black Ops Electric Purple Black Pearl

Peanut Butter Swirl Pink Flamingo Camo

Atlantis Global Gaze Rip Tide

○ = Circumference

△ = Width

□ = Distance

2 Diagonal circumference at apex of calcaneus

1 Smallest supra-malleolar circumference

3 Midfoot circumference at apex of arch

Width at met-heads (simulated weight bearing)

11 Posterior heel to finished lateral length

10 Posterior heel to finished medial length

LATERAL

MEDIAL

7


Posterior calcaneus to apex of 5th met-head

6

Posterior calcaneus to apex of 1st met-head

(Shoes are available for feet measuring 2 7/8" - 5 3/8")

4 Width at apex of malleoli



LATERAL

MEDIAL

9

Height to apex of lateral malleolus

8

Height to apex of medial malleolus

Wedging

Heel Neutral Medial _____" Lateral _____"

Forefoot Neutral Medial _____" Lateral _____"

Elevations

Heel Only Height _____

Heel & Sole Height _____

Special Modifications

Allard

Measurement required: Foot Length and Tibial Tuberosity and #10

KiddieGAIT Right Left Bilateral Soft Kit

KiddieROCKER Right Left Bilateral Soft Kit

Footwear

Shoes are available in sizes 3T to 12Y (measurement #6 = 2 7/8" - 5 3/8"). Please refer to our shoe form for volume orders.

PRODUCT SIZE COLOR

Athletic Shoes _____ White/Pink White/Blue Black

Sandals _____ Pink/Gray Black/Red

Toeless Socks _____ 2 pair/package

Shipping Options

UPS Ground No Charge

3 Day Select \$15.00

2nd Day (by 12:00 p.m.) \$25.00

Next Day Saver \$40.00

Next Day (by 10:30 a.m.) \$50.00

24 Hour Rush \$60.00 Added Fee

Shipping charges apply separately

Need by date: _____

Return Impressions Premodifications \$20.00/pr.

Send casts to: 17530 Dugdale Dr. South Bend, IN 46635