



# CFab KAFO Form

Version 01/2019

## Casts and Measurements Required

### ATTENTION:

### This is a 2 page form, please send both with order

Fax orders to 866-700-7837 | Customer Service 877-462-0711

17530 Dugdale Dr. South Bend, IN 46635 | www.surestep.net

#### PATIENT INFORMATION

Patient Last Name \_\_\_\_\_

Patient First Name \_\_\_\_\_

DOB \_\_\_\_\_ Weight \_\_\_\_\_ Sex M F

Dx \_\_\_\_\_

#### COMPANY INFORMATION

Company Name \_\_\_\_\_

Contact \_\_\_\_\_

Phone Number \_\_\_\_\_

Fax Number \_\_\_\_\_

Email \_\_\_\_\_

#### BILLING ADDRESS

City \_\_\_\_\_ ST \_\_\_\_\_ Zip \_\_\_\_\_

#### SHIPPING ADDRESS

City \_\_\_\_\_ ST \_\_\_\_\_ Zip \_\_\_\_\_

Order Confirmation:  Fax  Email

RIGHT  LEFT  BILATERAL

Date \_\_\_\_\_

Need by date \_\_\_\_\_

PO Number \_\_\_\_\_

#### SHIPPING OPTIONS

UPS Ground ..... No Charge

3 Day Select ..... \$25.00

2<sup>nd</sup> Day (by 10:30 a.m.) ..... \$30.00

Next Day Saver ..... \$45.00

Next Day (by 10:30 a.m.) ..... \$55.00

RUSH! ..... \$75.00 FEE

*SHIPPING CHARGES APPLY SEPARATELY*

#### NOTES

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

#### PATTERN OPTIONS

- Monsters
- Fairytale
- Peace Signs
- Comic Book
- Sports
- Butterflies
- Basketball
- Flamingos
- Pink Monkeys
- Spring Flowers
- Forest Friends
- Blossom
- Construction
- Dinosaurs
- Farm Friends
- Skateboard
- Police
- Pirates
- Nautical
- Busy Bees
- Traffic Soup
- Green Camo
- Pink Camo
- Happy Daisies
- Baseball
- Fairies
- Pink Leopard
- Music Notes
- Bugs
- White
- Yellow
- Pink
- Blue
- Cars
- Flames
- Monkeys
- Puppies

Other: \_\_\_\_\_

Non Corrosive Finish Color: \_\_\_\_\_

#### STRAP OPTIONS

- Black
- Yellow
- Green
- Royal Blue
- Purple
- Orange
- Lt. Blue
- White
- Pink
- Red

#### CHAFE OPTION

- Plastic (Standard)
- Dacron (Optional, no charge)



#### PLASTIC OPTIONS

- \_\_\_\_ Copolymer
- \_\_\_\_ Polypropylene
- \_\_\_\_ Polyethylene
- \_\_\_\_ Other
- Color \_\_\_\_\_
- Thickness \_\_\_\_\_

**FABRICATION**

- Plastic Double Upright KAFO** ..... X \_\_\_\_\_
- Plastic Single Upright KAFO** ..... X \_\_\_\_\_
- Plastic Static Knee KAFO** ..... X \_\_\_\_\_
- Material  Copoly  Polyeth  Polypro
- Thickness  1/8"  5/32"  3/16"  1/4"
- Pelvic Band (Sgl Up)** ..... X \_\_\_\_\_
- Pelvic Band (Dbl Up)** ..... X \_\_\_\_\_
- Quick-Release Proximal Section** ..... X \_\_\_\_\_
- Liner Above Knee** ..... X \_\_\_\_\_
- Thickness  1/8"  3/16"  1/4"
- Material  Volara
- Soft Pe-lite (White)
- P-Cell:  Pink  Black  Blue
- Puff Color \_\_\_\_\_
- Other \_\_\_\_\_
- Liner Below Knee** ..... X \_\_\_\_\_
- Thickness  1/8"  3/16"  1/4"
- Material  Volara
- Soft Pe-lite (White)
- P-Cell:  Pink  Black  Blue
- Puff Color \_\_\_\_\_
- Other \_\_\_\_\_
- Condyle Pad** ..... X \_\_\_\_\_
- Medial  Lateral
- Knee Pad (4 Buckle)** ..... X \_\_\_\_\_
- Knee Pad (5 Buckle)** ..... X \_\_\_\_\_
- Pre-Tibial Shell** ..... X \_\_\_\_\_
- Tuck-In  Overlap
- Thickness \_\_\_\_\_ Material \_\_\_\_\_
- Solid Ankle** ..... X \_\_\_\_\_
- Dual Adjustable Hinge** ..... X \_\_\_\_\_
- Surestep DA - Per Hinge, Per Device
- Molded Footplate, Stirrup Att.** ..... X \_\_\_\_\_
- Adjustable Plantarflexion Stop** ..... X \_\_\_\_\_
- SnapStop  Elite Stop
- Other \_\_\_\_\_
- T-Strap** ..... X \_\_\_\_\_
- Medial  Lateral
- Thickness \_\_\_\_\_
- Material  Duraflex  Polyethylene
- Molded Inner Boot** ..... X \_\_\_\_\_
- Inner Boot  SMO
- Thickness \_\_\_\_\_
- Material  Duraflex  Polyethylene  Volara
- Heel Post** ..... X \_\_\_\_\_
- Rigid Plastic  Crepe
- Dorsum Wings** ..... X \_\_\_\_\_
- Medial  Lateral

- MEDIAL LENGTH OF FOOT
- LATERAL LENGTH OF FOOT
- LENGTH TO 1ST MET HEAD
- LENGTH TO 5TH MET HEAD
- HIP CIR.
- HIP M/L
- WAIST CIR.
- WAIST M/L
- HEIGHT FROM TROC. TO WAIST



**COMPONENTS**

**Medial Knee Joint**

- MultiMotion  Ultraflex
- SS Free  HEKO
- Drop Lock  Ball Catch
- Other \_\_\_\_\_

**Lateral Knee Joint**

- MultiMotion  Ultraflex
- SS Free  HEKO
- Drop Lock  Ball Catch
- Other \_\_\_\_\_

**Medial Ankle Joint**

- MultiMotion  Ultraflex
- SS Free  COD
- DA Assist  Tam
- Other \_\_\_\_\_

**Lateral Ankle Joint**

- MultiMotion  Ultraflex
- SS Free  COD
- DA Assist  Tam
- Other \_\_\_\_\_

**COMPONENTS CONTINUED**

- Hip Joint Left \_\_\_\_\_
- Hip Joint Right \_\_\_\_\_
- Variable Abduct Hip Joint  x \_\_\_\_\_

**ADDITIONAL MEASUREMENTS**

- Sole to Fibula Head \_\_\_\_\_
- Sole to Knee Center \_\_\_\_\_
- Knee Center to Trochanter \_\_\_\_\_
- Knee Center to Perineum \_\_\_\_\_

**CAST MODIFICATIONS**

- Ankle:  Neutral / 90  as Casted
- Heel: \_\_\_\_\_
- Ankle: \_\_\_\_\_
- Forefoot: \_\_\_\_\_
- Knee:  Neutral / 180  as Casted
- Varus/Valgus: \_\_\_\_\_