

# STABILIZER



The Stabilizer is more than a preventative option for adults experiencing stability and balance issues. It is a solution that restores confidence, proprioception, and function. With the risk of falling greatly reduced, you're free to re-engage with day-to-day activities.

It features a low profile cuff as well as a molded inner boot. This prevents skin breakdown and keeps the heel back. And unlike other options, the Stabilizer allows greater sensory feedback by restricting the heel.

## INDICATIONS

Muscle weakness, ataxia, decreased balance response, gait abnormality, osteoarthritis (localized primarily in ankle/foot), arthropathy (ankle/foot), pain or instability of ankle/foot, mild drop foot, hemiplegia, Multiple Sclerosis

## CONTRAINDICATIONS

PTTD, charcot foot, spasticity

## NO CASTING REQUIRED

Significant time savings

4-day turnaround by measurements

## CONTACT US

17530 Dugdale Dr.

South Bend, IN 46635

P 877.462.0711 | F 866.700.7837

# STABILIZER

## PRODUCT FEATURES

Low profile cuff

Reinforcement for increased energy return

Molded Inner Boot provides consistent proprioceptive input through distribution of forces and circumferential compression

Shortened footplate for ease of donning in most any shoe



  
surestep

