



The **Stirrup AFO** is a custom ankle foot orthosis designed to treat chronic conditions of the foot and ankle, especially for patients needing more ankle stability and hindfoot control.

Utilizing industry-leading technology, the Surestep Stirrup AFO can be fabricated in 5 business days via cast or scan. The sleek, low-profile footplate and uprights are contoured to the patient's foot and ankle; helping to stabilize rotational forces at the midtarsal, subtalar, and talocrural joints.

The Stirrup AFO can be fabricated with either a solid or articulated ankle.

### Hinge Options

Ultra low profile | Surestep Free Motion  
3D Printed Polyurethane | Free Motion or Dorsi Assist



- » 3D printed product
- » 100% custom fabricated
- » 5-day turnaround
- » Company branding options

## CLINICAL INDICATIONS

- » Mild to moderate PTTD
- » Ankle, subtalar, or midtarsal instability
- » Mild ankle arthritis or DJD
- » Drop foot

### CONTACT US

17530 Dugdale Dr.  
South Bend, IN 46635  
P 877.462.0711  
F 866.700.7837

### Order Form

