



> GAUNTLET

The Surestep Gauntlet Series offers a superior approach to traditional custom leather ankle support systems. The combination of soft and rigid plastics allows for 3 unique versions to optimally treat severe (acute), moderate, or mild (early stage) conditions. Surestep Gauntlets are low profile and easy to fit.

INDICATIONS

Mild Gauntlet


Mild PTTD, mild ankle arthritis or DJD, ankle, subtalar, or midtarsal instability, mild sprains

Moderate Gauntlet

Moderate PTTD, moderate ankle arthritis or DJD

Severe Gauntlet

Severe PTTD, severe ankle arthritis or DJD



> STABILIZER

The Surestep Stabilizer is a solution that restores confidence, proprioception, and function. It features a low-profile cuff as well as a molded inner boot. This prevents skin breakdown and keeps the heel seated. And unlike other options, the Stabilizer allows greater sensory feedback through a unique open heel design. This allows the heel to 'feel' the ground increasing proprioceptive response.

INDICATIONS

Muscle weakness, ataxia, decreased balance response, gait abnormality, osteoarthritis (localized primarily in ankle/foot), arthropathy (ankle/foot), pain or instability of ankle/foot, mild drop foot, hemiplegia, MS



> STIRRUP

The Surestep Stirrup is a custom, ankle foot orthosis designed to treat chronic conditions of the foot and ankle, especially for patients needing more medial/lateral ankle stability and hindfoot control. The low-profile footplate and uprights are contoured to the patient helping to stabilize triplanar deformities.

INDICATIONS

Adult acquired flatfoot secondary to PTTD, muscle weakness of the ankle or foot, DJD of ankle or rearfoot, peroneal tendinopathy, drop foot



> FOOT ORTHOSIS

Surestep's 3D printed, functional foot orthoses are used to reduce abnormal motion and improve positioning for adult patients. With better foot and ankle alignment, patients will experience reduced pain, tenderness, and fatigue. The sulcus length top cover provides extra comfort and shock absorption.

INDICATIONS

Pes planus (flat feet), plantar fasciitis, heel spurs, degenerative joint disease, shin splints, posterior tibial tendonitis




> DCO

The Dynamic Cervical Orthosis is a solution offered for the treatment of head drop. It supports a neutral, midline position, while allowing normal movement. Resistance bands facilitate the strengthening of the musculature of the neck, while assisting with re-centering in midline. *Pediatric version available.*

INDICATIONS

Neck weakness, head drop (ptosis), Parkinson's disease, multiple sclerosis, myopathy, ALS, poor head control due to hypotonic cerebral palsy



adult
product line



Submitting orders is easy!

To submit an order, simply fill out a measurement form and send it to Surestep customer service representatives. Once your order is received, it will be processed and scheduled for fabrication.

Order forms are available on our website: surestep.net/clinicians
Email completed forms to orders@surestep.net or fax them to 866-700-7837

> Shipping

Surestep will ship via UPS ground, unless otherwise specified. For custom devices, ground shipping is complimentary. However, we do charge for upgraded shipping. Shipping charges apply for all off-the-shelf orders.

> Remake Policy

Surestep has a 90 day remake policy for all custom devices. A remake is warranted for any manufacturing defect or in the case of sudden growth or inadequate fit and function. Surestep reserves the right to charge additional fees when deemed necessary.

> Return Policy

Custom fabricated devices are non-returnable.

Our goal is to produce the best outcomes for each individual. If the original device ordered does not suit your patient's need, we will fabricate a different Surestep product that is more appropriate.

Any off-the-shelf items, including footwear, Derotation Straps, Prefab HEKOs, etc. may be returned within 14 days. Returns must be unworn, in the state you received them and in the original packaging. For detailed information on returning a specific product, please reference the instruction guide provided with your merchandise.



> Questions?

Surestep Customer Service Representatives

877.462.0711 | info@surestep.net

17530 Dugdale Dr. South Bend, IN 46635 | surestep.net

surestep

adult

product line

