

PDAC VERIFIED

The Stabilizer is more than a preventative option for adults experiencing stability and balance issues. It is a solution that restores confidence, proprioception, and function. With the risk of falling greatly reduced, you're free to re-engage with day-to-day activities.

It features a low profile cuff as well as a molded inner boot. This prevents skin breakdown and keeps the heel seated. And unlike other options, the Stabilizer allows greater sensory feedback through an open heel design.



INDICATIONS

Muscle weakness, ataxia, decreased balance response, gait abnormality, osteoarthritis (localized primarily in ankle/foot), arthropathy (ankle/foot), pain or instability of ankle/foot, mild drop foot, hemiplegia, Multiple Sclerosis

CONTRAINDICATIONS

PTTD, charcot foot, spasticity

NO CASTING REQUIRED

Significant time savings
4-day turnaround by measurements

■ **1940** Ankle foot orthosis, plastic or other material, custom-fabricated

■ **2280** Addition to lower extremity, molded inner boot

■ **2330** Addition to lower extremity, lacer molded to patient model, for custom fabricated orthosis only

■ **2820** Addition to lower extremity orthosis, soft interface for molded plastic, below knee section

CONTACT US

17530 Dugdale Dr.

South Bend, IN 46635

P 877.462.0711 | F 866.868.1890

STABILIZER

PRODUCT FEATURES



Low profile cuff

Molded Inner Boot provides consistent proprioceptive input through distribution of forces and circumferential compression

Shortened footplate for ease of donning in most any shoe




surestep

