



pediatric
custom product line




> **SMO/BIG SHOT SMO**

Surestep's dynamic SMO uses compression to stabilize the foot and ankle, allowing for the development of intrinsic muscles, muscle strategies, and movement patterns. The Big Shot SMO is designed with an inner boot for active children who still need the stability of an SMO. Cutouts allow for bony prominences.

INDICATIONS

Pronation, hypotonia, triplanar instability in weight bearing, inability to stand independently, mild toe walking, bony prominences (Big Shot only)




> **ADVANCED AFO**

Fixed ankle AFO easily transitions to a standard Surestep SMO.

INDICATIONS

Pronation, hypotonia, joint hypermobility, delayed standing, inability to stand for prolonged times



> **PULLOVER AFO**

Articulating dorsiflexion assist AFO that can be converted to a SMO by removing the proximal strut. Plantar flexion stop not available.

INDICATIONS

Pronation, hypotonia, mild-hemiplegia, drop foot, weakened dorsiflexors, sagittal plane instability




> **HINGED AFO**

An articulating AFO with free motion hinges or dorsiflexion assist hinges that has an adjustable plantar flexion stop.

INDICATIONS

Pronation, supination, triplanar instability, mild toe-walking, genu-recurvatum, mild foot drop



> **BIG SHOT AFO**

Fixed ankle AFO that provides relief for bony prominences.

INDICATIONS

Pronation, supination, hypotonia, hypertonia, limited weight bearing




> **INDY 2 STAGE**

Dynamic AFO/SMO that provides triplanar stability without restricting usage of intrinsic musculature of the foot.

INDICATIONS

Hypotonia, hypertonia, flexible pronation or supination, poor proprioceptive awareness, difficulty with transitional skills, sagittal and/or frontal plane weakness or instability, mild to moderate spasticity




> **HEKO**

Hyperextension knee orthosis with 4-axis hinge to stabilize the knee in all three planes while allowing full flexion and extension. Hinges can be sold separately.

INDICATIONS

Knee Instability, genu-recurvatum, flexible genu-valgum or genu-varum




> **TLSO**

Extremely flexible custom molded trunk stabilization orthosis.

INDICATIONS

Inability to control position in wheelchair, postural instabilities of the trunk, flexible kyphosis or lordosis





Submitting orders is easy!

To submit an order, simply fill out a measurement form and send it to Surestep customer service representatives. Once your order is received, it will be processed and scheduled for fabrication.

Order forms are available on our website: surestep.net/clinicians
Email completed forms to orders@surestep.net or fax them to **866.868.1890**

> Shipping

Surestep will ship via UPS ground, unless otherwise specified. For custom devices, ground shipping is complimentary. However, we do charge for upgraded shipping. Shipping charges apply for all off-the-shelf orders.

> Remake Policy

Surestep has a 90 day remake policy for all custom devices. A remake is warranted for any manufacturing defect or in the case of sudden growth or inadequate fit and function. Surestep reserves the right to charge additional fees when deemed necessary.

> Return Policy

Custom fabricated devices are non-returnable.

Our goal is to produce the best outcomes for each individual. If the original device ordered does not suit your patient's need, we will fabricate a different Surestep product that is more appropriate.

Any off-the-shelf items, including footwear, Derotation Straps, Prefab HEKOs, etc. may be returned within 14 days. Returns must be unworn, in the state you received them and in the original packaging. For detailed information on returning a specific product, please reference the instruction guide provided with your merchandise.



> Questions?

Surestep Customer Service Representatives

877.462.0711 | info@surestep.net
17530 Dugdale Dr. South Bend, IN 46635 | surestep.net



pediatric
custom product line

