


PATIENT

Last Name/Patient ID

First Name

Age Male Female

Weight Height

Dx

Right Left Bilateral

COMPANY

Name

Contact Name

Phone Fax

E-mail

BILLING / SHIPPING

BILLING

Address

City

ST/Prov Zip

Shipping address same as billing? Yes No

SHIPPING

Address

City

ST/Prov Zip

PO #

Order confirmation: Fax E-mail

FedEx SHIPPING OPTIONS

Need by Date

Ground®FREE

2 Day® A.M. 10:30 a.m. or noon 2nd business day\$30.00

Priority Overnight® 10:30 a.m. next business day\$55.00

First Overnight® 8:00 a.m. next business day\$100.00

Rush Fab 4 Day Turnaround \$100.00

CASTS

Sending Casts: Yes No

Return impressions pre-modifications \$20.00/pr.

Send casts:
 17530 Dugdale Drive
 South Bend, IN 46635

MEASUREMENTS

● Circumference ▲ Width ■ Distance

MEASURING VIDEO

CLICK OR SCAN

11 10

7 6

LATERAL MEDIAL

12 4

Circumference at distal fibula head

13

Desired finished height of AFO

9 8

* If number 5 measurement is 3" or greater, a cast is required.

DEVICE STYLE

SMO

Floor Reaction AFO

Solid Ankle AFO

Posterior Leaf AFO

Hinged AFO, Custom Fabricated

HINGE

Surestep Free Motion Surestep Dual Adjustable

Tamarack Free Motion Tamarack Dual Adjustable

Gaffney

Other:

STOPS

Surestop Snapstop

Elite Plastic

Other:

CORRECTION

Cast Modifications Neutral/90 As casted

Heel: Ankle: Forefoot:

PLASTIC OPTION

Polypropylene Polyethylene Copolymer Other:

Thickness: Color:

PATTERN CHAFE (D-RINGS)

Plastic Standard Dacron Optional, no charge

STRAP COLOR DORSAL PAD

Dorsal Pad Standard Neoprene Pad Optional, no charge

ADDITIONS/MODIFICATIONS (Additional charges may apply, see price list)

Liner

Volara (white) Puff

Other:

1/8" 3/16" 1/4"

Dorsum Wing

Medial Lateral

Heel Post

Rigid Plastic Crepe

Carbon Fiber Reinforcement

T-Strap

Medial Lateral

Duraflex Polyethylene

Condyle Extension

Medial Lateral

Molded Inner Boot

MPE Duraflex Foam

1/16" 3/32" 1/8"

TRIMLINE

Height: SMO Full

Boot: Short Surestep Full Length Footplate

Pre-Tibial Shell

Internal (tuck-in)

MPE Duraflex

3/32" 1/8"

External (overlap)
 Material same as device

Specify thickness:

Walking Base (Attached to AFO)

Open Heel

NOTES