



17530 Dugdale Dr. , South Bend, IN 46635

ARTICULATED AFO FORM 1970-ART-LE

Date / /

SUBMIT FORMS

Fax 866.868.1890 | E-mail orders@surestep.net

Customer Service 877.462.0711

PATTERNS
& STRAPS
CLICK OR SCAN



PATIENT

Patient ID		
Age	Weight	Height
Dx		
☐ Bilateral	☐ Left	☐ Right
Sending Casts: ☐ Yes ☐ No		

COMPANY

Name	
Contact Name	
Phone	Fax
E-mail	

SHIPPING

SHIPPING	Address	
	City	
	ST/Prov	Zip
PO #		
Order confirmation: ☐ Fax ☐ E-mail		

FedEx SHIPPING OPTIONS

Need by Date

☐ Ground®	FREE
☐ 2 Day® A.M. 10:30 a.m. or noon 2nd business day	\$30.00
☐ Priority Overnight® 10:30 a.m. next business day	\$55.00
☐ First Overnight® 8:00 a.m. next business day	\$100.00
Rush ☐ Rush Fabrication.....	\$100.00

FOOTWEAR	SIZE	QTY	COLOR
☐ Toddler Styles 3T - 9T			
☐ Youth Styles 10Y - 3Y			
☐ Heel Entry Style 3T - 3Y			White/Beige
☐ Sandals 3T - 3Y			Green

AFO SOCKS 1 PAIR/PACK	QTY	COLOR
☐ Infant 1 3/4" x 7 1/4"		
☐ Toddler 2" x 9"		
☐ Small 2 1/4" x 12"		
☐ Medium 3" x 15"		


MEASUREMENTS

● Circumference ▲ Width ■ Distance

MEASURING VIDEO



CLICK OR SCAN



* if number 5 measurement is 2 3/4" or greater, a cast is required.

The two standard AFO designs are listed below along with their standard modifications. If you would like to make changes, please select the base device and choose the options under the modifications columns that differ from the standards.

	☐ Functional Support	☐ Alignment Control	Modifications
Hinge	Surestep Free Motion	Tamarack Free Motion	☐ Surestep Free Motion
			☐ Tamarack Free Motion
Plastic	MATERIAL	Copolymer	Polypropylene
	THICKNESS	1/8"	1/8"
Stop	Plastic Stop	Plastic Stop	☐ Snapstop ☐ Surestop
			☐ Elite Plantarflexion Stop
Footplate Trimline	Full Length	Met Head	☐ Full Length ☐ Met Head
AFO Dorsum Wings	No	No	☐ Yes
Molded Inner Boot	MATERIAL	OP Flex	Low Density Polyethylene
	THICKNESS	1/8"	3/32"
	TRIM LINE (if applicable)	Full Length Footplate	Full Length Footplate
Ankle Pad	Yes	Yes	☐ No
Calf Padding	None	1/8" Puff	☐ None ☐ 1/8" Puff
Calf Strap	Layover Strap	Layover Strap	☐ Hook and Loop with Chafe
Chafe	Plastic	Plastic	☐ Dacron
Dorsal Padding	Dorsal Pad	Dorsal Pad	☐ Neoprene Pad

CORRECTION

Cast Modifications	☐ Neutral/90	☐ As casted
Heel:	Ankle:	Forefoot:

PATTERN	STRAP COLOR

ADDITIONS (Additional charges may apply, see price list)

☐ Non-Skid Sole	☐ Houdini Straps	☐ Adjustalift
☐ Metatarsal Pad	☐ Toe Plateau	

NOTES