



17530 Dugdale Dr. , South Bend, IN 46635

DASS NIGHT STRETCHING AFO

Dynamic Adjustable Stretching Strap

Date / /

SUBMIT FORMS

Fax 866.868.1890 | E-mail orders@surestep.net

Customer Service 877.462.0711

PATTERNS
& STRAPS
CLICK OR SCAN



PATIENT

Patient ID		
Age	Weight	Height
Dx		
☐ Bilateral	☐ Left	☐ Right
Sending Casts: ☐ Yes ☐ No		

COMPANY

Name		
Contact Name		
Phone	Fax	
E-mail		

SHIPPING

SHIPPING	Address		
	City		
	ST/Prov	Zip	
PO #			

Order confirmation: ☐ Fax ☐ E-mail

FedEx SHIPPING OPTIONS

Need by Date

☐ Ground		FREE
☐ 2 Day	A.M. 10:30 a.m. or noon 2nd business day	\$30.00
☐ Priority Overnight	10:30 a.m. next business day	\$55.00
☐ First Overnight	8:00 a.m. next business day	\$100.00

Rush ☐ Rush Fabrication.....\$100.00

FOOTWEAR	SIZE	QTY	COLOR
☐ Toddler Styles 3T - 9T			
☐ Youth Styles 10Y - 3Y			
☐ Heel Entry Style 3T - 3Y			White/Beige
☐ Sandals 3T - 3Y			Green

AFO SOCKS 1 PAIR/PACK

	QTY	COLOR
☐ Infant 1 3/4" x 7 1/4"		
☐ Toddler 2" x 9"		
☐ Small 2 1/4" x 12"		
☐ Medium 3" x 15"		


MEASUREMENTS

● Circumference ▲ Width ■ Distance

MEASURING VIDEO



CLICK OR SCAN



* if number 5 measurement is 2 3/4" or greater, a cast is required.

The Standard features are listed below Alignment Control. If you would like to specify modifications, select **Modifications** and customize. If option is left blank, the standard option will be used. Additional charges may apply, see price list.

	☐ Alignment Control	Modifications	
Medial Ankle Hinge	Goni-Joint	☐ Surestep Free Motion	
		☐ Tamarack	
Lateral Ankle Hinge	Pediatric Dual Adjustable Joint	☐ Surestep Free Motion	
		☐ Tamarack	
Plastic	MATERIAL	Copolymer	
	THICKNESS	3/32" (below 70lbs) ☐ 1/8" (over 70 lbs) ☐ Other:	
Molded Inner Boot	MATERIAL	OP Flex	
	THICKNESS	3/32"	
	TRIM LINE (if applicable)	Full Length Footplate	
Liner	MATERIAL	Volara (white)	
	THICKNESS	3/16"	
Walking Base	1/4" Soft PUFF midsole with thin tread		☐ Other:
Chafe (D-Rings)	Plastic		☐ Dacron (Optional, no charge)

CORRECTION

Cast Modifications	☐ Neutral/90	☐ As casted
Heel:	Ankle:	Forefoot:

PATTERN STRAP COLOR

NOTES	